

| ACCESSORIES |



Z.I. de Lambermont -151, rue des Ormes - B-4800 Verviers - BELGIUM
☎ : +32 (0)87 29 22 22 - 📠 : +32 (0)87 29 22 23
info@euromi.com - www.euromi.com

Your distributor

| CANNULAS |



INFILTRATION



Ø	LENGTHS	TYPE OF CANNULAS
2 mm ●	- / - / 200 mm / 250 mm / 300 mm	 Available in STAINLESS STEEL only
2.5 mm ●	- / - / 200 mm / 250 mm / 300 mm	
3 mm ●	- / - / - / 250 mm / 300 mm	

● The length in grey corresponds exclusively to custom-made cannula.

ASPIRATION


LIPOMATIC ASPIRATION

Available in STAINLESS STEEL or TUNGSTEN



CONVENTIONAL ASPIRATION


Available in STAINLESS STEEL or TUNGSTEN



TYPES OF CANNULAS

LIPOMATIC AND CONVENTIONAL ASPIRATION

Available in STAINLESS STEEL or TUNGSTEN




Ø	LENGTHS
3 mm ●	- / 150 mm / 200 mm / 250 mm / 300 mm / 350 mm
4 mm ●	- / 150 mm / 200 mm / 250 mm / 300 mm / 350 mm
5 mm ●	- / - / 200 mm / 250 mm / 300 mm / 350 mm
6 mm ●	- / - / - / 250 mm / 300 mm / 350 mm

● Lengths in grey correspond exclusively to custom-made cannulas

LIPOMATIC ASPIRATION ONLY

Only available in STAINLESS STEEL



BAYONET

Each type of cannula can be ordered custom-made in "bayonet" model




Ø	L
3 mm ●	250 mm
4 mm ●	300 mm





Any LIPOMATIC and CLASSICAL aspiration cannulas are available in TITANIUM upon request only.

LIPOFILLING



Ø	LENGTHS	TYPES OF CANNULAS
2 mm ●	100 mm / 150 mm / 200 mm	 1 B
2.5 mm ●	100 mm / 150 mm / 200 mm	
3 mm ●	100 mm / 150 mm / 200 mm	
2.5 mm ●	100 mm / 150 mm / 200 mm	 MERCEDES
3 mm ●	100 mm / 150 mm / 200 mm	
3 mm ●	150 mm / 200 mm	
		 VITERBO

Ø	LENGTHS	TYPES OF CANNULAS
2.5 mm ●	150 mm / 200 mm	 ABBOUD I
3 mm ●	150 mm / 200 mm	
3 mm ●	150 mm / 200 mm	 ABBOUD P

INFILTRATION

ASPIRATION

LIPOFILLING